











	AA battery-option			Rechargeable battery-option						
	Ameda Mya Joy	Lansinoh SmartPump 2.0	Medela PIS with MaxFlow ¹	Cimilre S6+	Cimilre S5+	Cimilre P1	Medela Freestyle Flex	Medela Swing Maxi	Spectra S1 Plus	Zomee Z1
Breast Pump Features										
Item ID	131A40	53216	101041360	S6P-100-00-00	S5P000	P1-100-00-00	101037980	101043614	SP08-L1	ZOMEE Z1
Weight	0.5 lbs	3 lbs	1.2 lbs	1.9 lbs	0.9 lbs	0.55 lbs	~1 lbs	<1 lbs	4 lbs	3 lbs
Flange Sizes	25mm, 28.5mm	25mm	24mm	24mm, 28mm	24mm, 28mm	24mm, 28mm	21mm, 24mm	21mm, 24mm	24mm, 28mm	24mm, 27mm
Warranty	2 yr	1 yr	1 yr	2 yr	2 yr	2 yr	1 yr	1 yr	2 yr	1 yr
Max Suction Strength²	-250 mmHg	-250 mmHg	-295 mmHg	-285 mmHg	-295 mmHg	-285 mmHg	-270 mmHg	-270 mmHg	-270 mmHg	-260 mmHg
Skinny Neck Bottles	✓		✓				✓	✓		
Wide Neck Bottles		✓		✓	✓	✓			✓	✓
Bluetooth		✓					✓	✓		
Digital Display	✓	✓		✓	✓	✓	✓		✓	✓
Personalize				✓	✓	✓				
Noise volume²	<45 dB	Not available	Not available	<45 dB	<45 dB	<45 dB	Not available	Not available	<45 dB	≤40 dB
Night-Light				✓					✓	
Dual Motor					✓					
Independently Adjustable Speed & Suction				✓					✓ ³	
Letdown/Expression Mode(s)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

¹ External battery pack. ² Provided by manufacturers. ³ Only in Expression Mode.

	Wearable			Stationary	
	Elvie Double Electric Breast Pump	Elvie Stride Hands-Free Breast Pump	Willow [®] Generation 3 Wearable Double Hands-Free Electric Breast Pump	Cimilre S6	Spectra S2 Plus

Breast Pump Features



Item ID	EP01-02	EB01-02	PDW48	S6-100-00-00	SP06-L1
Weight	~1lb combined	~1lb combined	~1lb combined	1.7 lbs	4 lbs
Flange Sizes	24mm, 28mm	24mm	24mm	24mm, 28mm	24mm, 28mm
Warranty	2 yr	2 yr	1 yr	2 yr	2 yr
Max Suction Strength²	-220 mmHg	-300 mmHg	-245 mmHg	-285 mmHg	-270 mmHg
Skinny Neck Bottles	✓				
Wide Neck Bottles				✓	✓
Bluetooth	✓	✓	✓		
Digital Display	✓ ⁴	✓ ⁴	✓ ⁴	✓	✓
Personalize				✓	
Noise volume²	No discernable dB level	No discernable dB level	Not available	<45 dB	≤45 dB
Night-Light				✓	✓
Dual Motor					
Independently Adjustable Speed & Suction				✓	✓ ³
Letdown/Expression Mode(s)	✓	✓	✓	✓	✓

² Provided by manufacturers.

³ Only in Expression Mode. ⁴ Digital display is on mobile app.